

Acorn 120 Superglide



BABYLON SURGICAL SUPPLIES, INC. • DEER PARK, NY 11729 • (631) 242-7900 • www.babylonsurgicalsupplies.com



Swivel seat puts users feet firmly on to the top landing area making getting off the seat easier and safer



Acorn 120 Superglide. Remote controls allow users to 'call' or 'send' the lift up and down the staircase. We provide two remote control units with each lift.



Acorn 120 Superglide. Our focus on safety features ensures we provide the best seatbelts for our customers.



Acorn 120 Superglide. Multiple safety edges ensure that the stairlift will stop immediately if an obstruction is encountered on the staircase.



Acorn 120 Superglide. A lockable on/off switch allows the user to prevent others (e.g. children) from using the lift. It also has a digital diagnostic display.



Acorn 120 Superglide.

Track fits directly on to stair nose and fixes on to stair tread



Acorn 120 Superglide. The slimline fold-away design means that other members of the family can still use the stairs as normal.



Acorn 120 Superglide. The directional paddle switches are operated by the lightest touch, even by those with limited dexterity.




Brooks



The stairlift fixes quickly and directly on to your staircase, not to the wall.



Technical Specification for 120 Superglide



DIMENSIONS

	mm	inch
A - Width outside armrests	560	22
B - Width inside armrests	445	17.5
C - Width of backrest	370	14.5
D - Depth of seat	410	16
E - Top of footrest to top of seat	445	17.5
F - Required space from bottom of stair	660	26
G - Wall to Seat width	40	1.5
H - Overall height of seat	495	19.5
I - Overall height	1000	40
J - Footrest surface to floor	115	4.5
K - Overall depth (folded)	325	13
L - Overall depth (footrest down)	585	23

WIDE SEAT OPTION

B - Width inside armrests	520	20
F - Required space from bottom of stair	710	28
D - Depth of seat	410	16

SPECIFICATIONS

Motor Power	0.5 hp
Travel Speed	No greater than 0.15m/s
Method of Drive	Rack and Pinion
Motor Output Torque	108 Nm
Power Supply	24v DC (battery)
Maximum Capacity	294 lbs (optional 350lbs)
Track	Extruded Aluminium
Mains Supply	110v, 1.5amp Charger

ACORN, in keeping with its policy of continual development, reserves the right to change specifications without notice. All measurements are approximate.



Perch Seat and Sit/Stand stairlifts



Acorn 120 Perch Stairlift.
Ideal for when the user
may have trouble sitting
or bending at the knee
and therefore cannot use a
standard stairlift.



Acorn 120 Sit/Stand Stairlift.
For multiple users with
different requirements.
This model allows the user
to choose whether to sit
or stand.